


EASY TO USE EASY TO TRUST



ASFV qPCR Diagnostic Kit

Multiplex solution
in diagnostics



Introducing ASFV qPCR Diagnostic Kit, the easy solution for everyday veterinary diagnostics.

In the field of diagnostics, **accuracy and reliability** are essential. The TargetEx ASFV qPCR Diagnostic Kit empowers laboratories with fast and confident detection of African Swine Fever Virus thanks to its **advanced triplex qPCR technology**.

Delivering **results in 45-50 minutes**, the kit facilitates efficient testing, reduces uncertainty, and enables **confident decision-making** when time is of the essence.

Independently **validated** by the EU Reference Laboratory for ASFV, this assay ensures **reliable and consistent** performance.

Designed for **maximum flexibility**, the assay delivers robust performance across a **wide variety of sample types**.

FEATURES AT A GLANCE:

Triplex qPCR: High-performance technology with all relevant controls included in a single kit (FAM, HEX and Cy5 channels).

Fast Results: Rapid protocol delivers reliable results in just 45 minutes.

Quality You Can Trust: High analytical and diagnostic specificity, low limit of detection, and certified under ISO 9001:2015, with validation by the INIA-CSIC EURL for ASFV.

Wide sample compatibility: Works perfectly with serum, whole blood, and various organ and tissue samples (e.g. bone marrow), with the option for sample pooling.

Versatility: Tested on a wide range of qPCR thermocyclers, e.g. BioRad CFX96, Rotor-Gene Q, Agilent AriaMX, Quant Studio 5.

Ready to Simplify Your Diagnostic Procedures?



Contact us for a sample & unlock the power of ASFV qPCR Kit today.

ABOUT US

TargetEx is the powerful ally
in your **biotechnological**
advancements

TargetEx is your trusted ally in Biotechnology Research and Development, dedicated to discovering the best approach for your unique challenges. Our mission is to build collaborative partnerships and find bespoke solutions, ensuring that our valued clients not only find support but also build long-lasting relationships.



Let's start
a **conversation!**

M: info@targetex.com

T: +36 30 252 0971



Follow us on LinkedIn



targetex.com